

Cleanroom Pharmacy Syringe Packs

- › Micronclean Designed: Specifically for NHS compounding use
- › Radiation stable European polymer
- › Highest quality graduations in millilitres only
- › Compatible with most syringe drivers
- › Luer lock, complete with tip protector
- › Offers significant cost savings


» Ensuring Compliance

Compliance 100[®] products are supplied in transparent, easy-tear, quad laminate pouches which are VHP, alcohol and peracetic acid impermeable, rendering them suitable for all transfer processes

Overall Benefits Table	
Feature	Benefit
Multiple syringe packs	Increased efficiency of transfer times and batch paperwork, leaving more technician time to focus on aseptic production itself.
Packaging is all transparent	Rapid intuitive confirmation of contents and QC checking prior to transfer. All packaging is non-shedding.
Quad laminate pouch	Transparent, strong, easy to tear, with burst resistant broad seal and rounded corners to avoid puncturing outer bag.
Pouch size optimal	Pouch size is designed to meet the dual requirements of coping with the volume of syringes whilst avoiding spillage of syringes on opening in the LAF or Isolator.
Single batches of syringes	All pouches contain only one batch no. of syringes, streamlining and simplifying your batch paperwork to a minimum.
100% tipped syringes	The luer lock syringe tip protector ensures the sterility of the vulnerable tip during manipulation.
Consolidated Certificate of Conformity	A single document which is easy to file and cross reference.
Large stocks held by Micronclean	Security of supply.
Logistics by Micronclean	Easy to set up service with weekly free deliveries.


Key Features & Benefits

- › Sterile tip protector in unique colour
- › Latex free bung
- › Extremely tight siliconisation controls
- › Crystal clear indelible graduations in millilitres only
- › European gamma irradiation stable polymer
- › Strong 'anti-pull out' mechanism
- › Retrofits into most popular syringe drivers



» Pouched Product

All product is provided with a second polyethylene over-pouch bag fitted



» Easy tear opening

» Clear VHP, alcohol and peracetic acid impermeable pouch




Ordering Information

Micronclean Syringe Size

	3ml	5ml	10ml	20ml	30ml	50ml
Pouch of 5 order code	03ML05MC	05ML05MC	10ML05MC	20ML05MC	30ML05MC	50ML05MC
No. of pouches per box	20	20	20	20	20	20
Pouch of 10 order code	03ML10MC	05ML10MC	10ML10MC	20ML10MC	30ML10MC	50ML10MC
No. of pouches per box	10	10	10	10	10	10
Pouch of 25 order code	03ML25MC	05ML25MC	10ML25MC	20ML25MC	30ML25MC	50ML25MC
No. of pouches per box	10	10	10	10	10	10

Right: Clear graduations.

Far right: Unique colour tip protector.





Easy to Deal With

Contact Micronclean

For any enquiries please contact Micronclean - we are happy to talk through any needs you may have.



Please contact us for a copy of our Single Use Consumables brochure.



Micronclean
Roman Bank
Skegness
Lincolnshire PE25 1SQ

Tel: 01754 767377
Fax: 0845 2992166

Email: enquiries@micronclean.co.uk

Web: www.micronclean.co.uk

Garments & Mops website

Web: www.micronclean.eu

Consumables website