

Engineered for Protection. Designed for Comfort.

KIMTECH PURE* G3 Sterile STERLING* Nitrile Gloves

KIMBERLY-CLARK* has developed a low-cost sterile nitrile glove combining the sensitivity of latex with the protection of nitrile. Our innovative STERLING* gloves manufacturing technology enables us to deliver a comfortable and affordable glove with all the protection and cleanliness our customers have come to expect.

When double-donned, the KIMTECH PURE* G3 Sterile STERLING* Nitrile Gloves deliver the equivalent feel & dexterity of a single latex glove without the risk of Type I latex allergic reactions. Protection has never been so comfortable.

The KIMTECH PURE* G3 Sterile STERLING* Nitrile Glove is also environmentally responsible. We have improved packaging to deliver 50% more gloves per case than traditional gloves. Average use of 25 cases per month results in reduction of 1 ton of waste per year vs. traditional latex gloves.

Calculate Your Saving:
www.kimtech.com/sterlingcalculator



KIMTECH
PURE* BRAND



Quality Standards


- PPE Category III product classified by EC Council Directive 89/686/EED
- Tested in accordance with EN Norms EN 420-2003
- Cleanroom packaged to meet the needs of Class 100 Cleanrooms
- Meets or exceeds AQL level of 1.5 for pinholes
- Manufactured in accordance with quality System ISO 9001

Cleanliness Characteristics

Parameter	Limit	Test Method
Particles		
Per cm ² ≥ 0.5 micron	1200	IEST-RP-CC005

Physical Properties

Characteristics	Value								Test Method
Freedom from Holes	1.5 AQL¹								ASTM D 5151
Tensile Properties	Tensile Strength		Ultimate Elongation						ASTM D 412 and ASTM D 573
Before Aging	42 MPa, nominal		650%, nominal						
After Accelerated Aging	38 MPa, nominal		550%, nominal						
Dimensional	Measured Point		mm					mil	ASTM D 3767 and ASTM D 6319
Nominal Thickness	Middle Finger		0.10					3.9	
Palm			0.08					3.1	
Cuff			0.07					2.8	
Nominal Length			305mm / 12 inches						
Palm Widths	6	6.5	7	7.5	8	8.5	9	10	ASTM D 3767 and ASTM D 6319
Nominal Width (mm)	80	87	94	98	109	114	120	128	
Endotoxin	(Maximum)								LAL Kinetic Turbidimetric Method
Endotoxin Units/pair	20								



KIMTECH PURE* G3 Sterile STERLING* Nitrile Gloves

Packaged in a Class 100 Cleanroom

Description	6	6.5	7	7.5	8	8.5	9	10	Pairs/Bags	Total/Case
KIMTECH PURE* G3 Sterile STERLING* Nitrile Gloves - 12"	11821	11822	11823	11824	11825	11826	11827	11828	30	300

Our Guarantee

Your total satisfaction means everything to us. If, for any reason, our products do not meet your expectations, Kimberly-Clark will reimburse you† for your initial purchase, via FREE product, for up to \$1,000. For more information on KIMBERLY-CLARK PROFESSIONAL*, visit us online at www.kcprofessional.com, contact your Kimberly-Clark Sales Representative, or call us at 1-888-346-GOKC (4652).

† Guarantee extended to consuming end-user accounts only.

Reduce Today, Respect Tomorrow is the KIMBERLY-CLARK PROFESSIONAL® approach to sustainability. It begins with the understanding that the way we use resources today shapes the world of tomorrow. And it has led us to focus on reducing consumption at every stage of the product lifecycle – from design and manufacture to distribution and disposal. Reduction is the key to lowering the environmental impact of our activities as well as those of customers. It has also been crucial in achieving our number one position five years in a row on the Dow Jones Sustainability Index. To learn more about Reduce Today, Respect Tomorrow and how we can reduce consumption in your business, visit www.kcreducetoday.com



Kimberly-Clark warrants that its products comply with K-C's standard specifications as of the delivery date to K-C's authorized distributors/direct purchasers. Except to the extent prohibited by applicable law: 1) **THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE**; 2) K-C is not liable for any kind of special, incidental, or consequential damages; and 3) K-C's liability for breach of contract, tort, or other cause of action shall not exceed the product purchase price. Purchasers and users are deemed to have accepted the above warranty and limitation of liability, and cannot change the terms by verbal agreement or by any writing not signed by K-C. To the extent required by applicable law, K-C does not limit its liability for death/injury resulting from K-C's negligence.

®/* Trademarks of Kimberly-Clark Worldwide, Inc. or its affiliates. Marques déposées de Kimberly-Clark Worldwide, Inc. ou de ses filiales. © 2009 KCWW. K01722 K3619-09-01 11/09

Effective November 2009