

**MICRONCLEAN COMPLIANCE 100™  
WIPES SPECIFICATION**

**Product Code:** MCZSLA50S

**Product Description:** Sterile cleanroom wipes made from single layer continuous filament polyester with laser cut edges impregnated with 70% isopropyl alcohol / 30% water for injection

**Pack Size:** 50 wipes / pack & 5 packs per case

**Pouch Format:** Micro25 and treble bagging

**Technical Specification**

**Fabric:** 100% polyester; single layer construction; 132 g/m<sup>2</sup> weight; tested to ASTM-D-629

**Wipe Size:** 300mm x 300mm

**Alcohol:** 70% USP IPA / 30% USP WFI

**Decontamination:** Aqueous washing incorporating thermal disinfection to minimum 71°C / 3 minutes; product packed in ISO 14644-1 Class 4 (at rest), EU GMP Grade B cleanroom

**Sterilisation:** Validated sterilisation by gamma irradiation (minimum 25 kGy)

**Shelf-Life:** Two years from manufacturing date

**Quality Assurance:** Product manufactured under BS EN ISO 9001:2000 registered systems