

Sterile Micronclean Pharmawipe (Dry) 12'' x 12''

Description

Dry polycellulose flat wipes packed in packs of 25 x 4 with 17 of these per pack (1700 wipes in grand total)

Cleanroom classification

Non-woven substrates are manufactured and packed in a controlled environment of Class 10,000 (ISO 7)

Sterilisation

Product is sterilised by gamma irradiation, min 25 kG max 45 kG in a process controlled by Micronclean.

Packaging Material

First bag = Polyethylene Second bag = Polyethylene