

Sterile Micronclean Pharmawipe (Dry) 9” x 9”

Description

Dry polycellulose c-fold wipes packed in bags of 10 x 10 with 27 of these per Box (2700 wipes in grand total)

Cleanroom classification

Non-woven substrates are manufactured and packed in a controlled environment of Class 10,000 (ISO 7)

Sterilisation

Product is sterilised by gamma irradiation, min25 kG max 45 kG in a process controlled by Micronclean.

Packaging Material

First bag = Polyethylene Second bag = Polyethylene